

Anti-Nivolumab Polyclonal Antibody

Catalog No: #V6869



Package Size: #V6869-1 50ug #V6869-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Nivolumab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Nivolumab
Conjugates	Unconjugated
Target Name	Nivolumab, BMS-936558, MDX-1106, MDX1106, MDX 1106, ONO-4538, OPDIVO
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.