

Anti-Atezolizumab Polyclonal Antibody

Catalog No: #V6872



Package Size: #V6872-1 50ug #V6872-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Atezolizumab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Atezolizumab
Conjugates	Unconjugated
Target Name	Rabbit polyclonal to Tecentriq.
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.