

# Anti-Ocrelizumab Polyclonal Antibody

Catalog No: #V6877



Package Size: #V6877-1 50ug #V6877-2 100ug


Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                     |
|-----------------------|-----------------------------------------------------|
| Product Name          | Anti-Ocrelizumab Polyclonal Antibody                |
| Host Species          | Rabbit                                              |
| Clonality             | Polyclonal                                          |
| Isotype               | IgG                                                 |
| Purification          | Purified by antigen affinity column.                |
| Applications          | ELISA                                               |
| Immunogen Description | Ocrelizumab                                         |
| Conjugates            | Unconjugated                                        |
| Target Name           | Rabbit polyclonal to Ocrevus.                       |
| Formulation           | 0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300. |

## Images

Anti-Ocrelizumab Non-Neutralizing Antibody (AD2133) binds with Ocrelizumab



Detects Ocrelizumab in indirect ELISA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.