

# Anti-Ofatumumab Polyclonal Antibody

Catalog No: #V6880



Package Size: #V6880-1 50ug #V6880-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                     |
|-----------------------|-----------------------------------------------------|
| Product Name          | Anti-Ofatumumab Polyclonal Antibody                 |
| Host Species          | Rabbit                                              |
| Clonality             | Polyclonal                                          |
| Isotype               | IgG                                                 |
| Purification          | Purified by antigen affinity column.                |
| Applications          | ELISA                                               |
| Immunogen Description | Ofatumumab                                          |
| Conjugates            | Unconjugated                                        |
| Target Name           | Rabbit polyclonal to HuMax-CD20.                    |
| Formulation           | 0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.