

Anti-Imdevimab Polyclonal Antibody

Catalog No: #V6882



Package Size: #V6882-1 50ug #V6882-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Imdevimab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Imdevimab
Conjugates	Unconjugated
Target Name	Rabbit polyclonal to REGN-10987.
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.