

Anti-Ramucirumab Polyclonal Antibody

Catalog No: #V6884

Package Size: #V6884-1 50ug #V6884-2 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Anti-Ramucirumab Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA
Immunogen Description	Ramucirumab
Conjugates	Unconjugated
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.