

Anti-Mugwort Art v 1 Polyclonal Antibody

Catalog No: #V7027



Package Size: #V7027-1 50ug #V7027-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Mugwort Art v 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Artemisia vulgaris (Mugwort)
Immunogen Description	E. coli - derived recombinant Mugwort Art v 1 (Ala25-His132).
Conjugates	Unconjugated
Accession No.	Q84ZX5
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Images



Western blot analysis was performed using anti-ARTV1 polyclonal antibody at 1ug/mL on various samples.Lane 1: recombinant Artemisia vulgaris Artv1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.