

Anti-Neosartorya fumigata Asp f 1/mitF Polyclonal Antibody

Catalog No: #V7028

Package Size: #V7028-1 50ug #V7028-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Neosartorya fumigata Asp f 1/mitF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Immunogen Description	E. coli - derived recombinant Neosartorya fumigata Asp f 1/mitF (Ala28-His176).
Conjugates	Unconjugated
Accession No.	P67875
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Images



Western blot analysis was performed using anti-ASP f1 polyclonal antibody at 1ug/mL on various samples.Lane 1: recombinant Aspergillus fumigatus Asp f 1/mitF

Note: This product is for in vitro research use only and is not intended for use in humans or animals.