

Anti-Neosartorya fumigata Asp f 15 Polyclonal Antibody

Catalog No: #V7029

Package Size: #V7029-1 50ug #V7029-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-Neosartorya fumigata Asp f 15 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (<i>Aspergillus fumigatus</i>)
Immunogen Description	E. coli - derived recombinant Neosartorya fumigata Asp f 15 (Ser34-Val152).
Conjugates	Unconjugated
Accession No.	O60022
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.