

Anti-European white birch Bet v 1/Allergen Bet v I-A Polyclonal Antibody

Catalog No: #V7030

Package Size: #V7030-1 50ug #V7030-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-European white birch Bet v 1/Allergen Bet v I-A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Betula pendula (European white birch)
Immunogen Description	E. coli - derived recombinant European white birch Bet v 1/Allergen Bet v I-A (Met1-Asn160).
Conjugates	Unconjugated
Accession No.	P15494
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.