

Anti-German cockroach Bla g 1/Allergen Bla g I Polyclonal Antibody

Catalog No: #V7033

Package Size: #V7033-1 50ug #V7033-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-German cockroach Bla g 1/Allergen Bla g I Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Blattella germanica (German cockroach)
Immunogen Description	E. coli - derived recombinant German cockroach Bla g 1/Allergen Bla g I (Arg34-Ser225).
Conjugates	Unconjugated
Accession No.	Q9UAM5
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.