

Anti-Blomia tropicalis (Mite) Blo t 5/BLOT5 Polyclonal Antibody

Catalog No: #V7036

Package Size: #V7036-1 50ug #V7036-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-Blomia tropicalis (Mite) Blo t 5/BLOT5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Blomia tropicalis (Mite)
Immunogen Description	E. coli - derived recombinant Blomia tropicalis (Mite) Blo t 5/BLOT5 (Gln18-Gln134).
Conjugates	Unconjugated
Accession No.	O96870
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.