

Anti-Bovine ALB Polyclonal Antibody

Catalog No: #V7041



Package Size: #V7041-1 50ug #V7041-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Bovine ALB Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Bovine
Immunogen Description	E. coli - derived recombinant Bovine ALB (Ser18-Pro327).
Conjugates	Unconjugated
Accession No.	P02769
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Images



Western blot analysis was performed using anti-ALB polyclonal antibody at 1ug/mL on various samples.Lane 1: Hela cell lysateLane 2: LO2 cell lysateLane 3: NIH3T3 cell lysateLane 4: Rat kidney lysateLane 5: Rat spleen lysate

Note: This product is for in vitro research use only and is not intended for use in humans or animals.