

Anti-European house dust mite Der p 10/Tropomyosin Polyclonal Antibody

Catalog No: #V7049

Package Size: #V7049-1 50ug #V7049-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-European house dust mite Der p 10/Tropomyosin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Dermatophagoides pteronyssinus (European house dust mite)
Immunogen Description	E. coli - derived recombinant European house dust mite Der p 10/Tropomyosin (Met1-Tyr284).
Conjugates	Unconjugated
Accession No.	O18416
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.