

Anti-Cat Fel d 1A/CH1/Allergen Fel dI Polyclonal Antibody

Catalog No: #V7054

Package Size: #V7054-1 50ug #V7054-2 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anti-Cat Fel d 1A/CH1/Allergen Fel dI Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Felis catus (Cat) (Felis silvestris catus)
Immunogen Description	E. coli - derived recombinant Cat Fel d 1A/CH1/Allergen Fel dI (Glu23-Cys92).
Conjugates	Unconjugated
Accession No.	P30438
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.