

Anti-Cat Fel d 1B/CH2/Allergen FdI Polyclonal Antibody

Catalog No: #V7055

Package Size: #V7055-1 50ug #V7055-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Cat Fel d 1B/CH2/Allergen FdI Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by antigen affinity column.
Applications	ELISA, IHC, WB
Species Reactivity	Felis catus (Cat) (Felis silvestris catus)
Immunogen Description	E. coli - derived recombinant Cat Fel d 1B/CH2/Allergen FdI (Val18-Arg109).
Conjugates	Unconjugated
Accession No.	P30440
Formulation	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

ELISA:1:5000-1:20000;IHC:1:50-1:500;WB:1:500-1:2000

Images



Western blot analysis was performed using anti-fel d 1 chain 2 polyclonal antibody at 1ug/mL on various samples.Lane 1: recombinant Cat Fel d 1B/CH2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.