

Anti-Ganglioside GM1 hIgM Recombinant Antibody

Catalog No: #V8060



Package Size: #V8060-1 50ug #V8060-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Ganglioside GM1 hIgM Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	GL2049
Isotype	IgM
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Human
Conjugates	Unconjugated
Target Name	Ganglioside GM1, Ganglioside A2, Ganglioside G4, Ganglioside GGtet1, Ganglioside GI, Monosialoganglioside GM1
Accession No.	ChEBI: 61048
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Ganglioside GM1 hIgM Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.