

Anti-Human MUC-1 (Tn/STn) Recombinant Antibody, APC

Catalog No: #V8068

Package Size: #V8068-1 50T #V8068-2 100T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Anti-Human MUC-1 (Tn/STn) Recombinant Antibody, APC
Host Species	Human
Clonality	Monoclonal
Clone No.	GL2089
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, FCM, WB
Species Reactivity	Human
Conjugates	APC
Target Name	Breast carcinoma-associated antigen DF3, CA 15-3, CD227, Cancer antigen 15-3, Carcinoma-associated mucin, EMA, Episialin, H23AG, KL-6, Krebs von den Lungen-6, MUC-1, MUC1, MUC1-CT, MUC1-NT, MUC1-alpha, MUC1-beta, Mucin-1, Mucin-1 subunit alpha, Mucin-1 subunit beta, PEM, PEMT, PUM, Peanut-reactive urinary mucin, Polymorphic epithelial mucin, Tumor-associated epithelial membrane antigen, Tumor-associated mucin
Accession No.	P15941
Formulation	0.01M PBS, pH 7.4, 50% glycerol, 0.05% Proclin 300.
Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.

Application Details

Purity:>90% as determined by SDS-PAGE.

Images



Flow-cytometry using APC anti-human CD227 antibody. MCF-7 cells were stained with an irrelevant antibody (Blue Histogram) or an APC anti-human CD227 monoclonal antibody (Yellow Histogram) at a concentration of 5 µg/ml for 30 mins at RT. After washing, and cells analysed on a NovoCyte Flow Cytometer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.