

Anti-Sahara scorpion AaH I/Neurotoxin I Recombinant Antibody

Catalog No: #V9003



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9003-1 50ug #V9003-2 100ug


Description

Product Name	Anti-Sahara scorpion AaH I/Neurotoxin I Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0024
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Androctonus australis (Sahara scorpion)
Conjugates	Unconjugated
Target Name	Alpha-mammal toxin AaH1, AaHI, Neurotoxin 1, Neurotoxin I
Accession No.	P01479
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Sahara scorpion AaH I/Neurotoxin I Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.