

Anti-Bacillus anthracis pagA/Protective antigen/PAD4
Recombinant Nanobody

Catalog No: #V9013



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9013-1 50ug #V9013-2 100ug


Description

Product Name	Anti-Bacillus anthracis pagA/Protective antigen/PAD4 Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	ST0021
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column.
Applications	ELISA
Species Reactivity	Bacillus anthracis
Conjugates	Unconjugated
Target Name	pagA, PA, Protective antigen, Anthrax toxins translocating protein, PA-83, PA83, Bacillus anthracis
Accession No.	P13423
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Bacillus anthracis PAD4 Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.