

Anti-Bacillus anthracis lef/LF/Lethal factor
Recombinant Antibody

Catalog No: #V9016



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9016-1 50ug #V9016-2 100ug


Description

Product Name	Anti-Bacillus anthracis lef/LF/Lethal factor Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0039
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, IP, SPR
Species Reactivity	Bacillus anthracis
Conjugates	Unconjugated
Target Name	lef, LF, Lethal factor, Anthrax lethal toxin endopeptidase component, Bacillus anthracis
Accession No.	P15917
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Bacillus anthracis lef/LF/Lethal factor Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.