

Anti-Bacillus anthracis cya/EF/Edema factor
Recombinant Antibody

Catalog No: #V9018



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9018-1 50ug #V9018-2 100ug


Description

Product Name	Anti-Bacillus anthracis cya/EF/Edema factor Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0043
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, SPR
Species Reactivity	Bacillus anthracis
Conjugates	Unconjugated
Target Name	Calmodulin-sensitive adenylate cyclase, 4.6.1.1, ATP pyrophosphate-lyase, Adenylyl cyclase, Anthrax edema toxin adenylate cyclase component, Edema factor, EF, cya, pXO1-122, BXA0141, GBAA_pXO1_0142, Bacillus anthracis
Accession No.	P40136
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Bacillus anthracis cya/EF/Edema factor Recombinant Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.