

Anti-Mexican scorpion Beta-toxin Cl13/Cb3 Recombinant Antibody

Catalog No: #V9022

Package Size: #V9022-1 50ug #V9022-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Mexican scorpion Beta-toxin Cl13/Cb3 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0023
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization, SPR
Species Reactivity	Centruroides (Mexican scorpion)
Conjugates	Unconjugated
Target Name	Beta-toxin Cl13, CliNaTxBet37, Beta-toxin Cb3, Centruroides limpidus (Mexican scorpion), Centruroides baergi (Scorpion) (Centruroides nigroviratus baergi)
Accession No.	C0HK69, C0HLR5
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Mexican scorpion Beta-toxin Cl13/Cb3 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.