

Anti-Clostridium botulinum HA-33/Hemagglutinin 33 Recombinant Antibody

Catalog No: #V9027

Package Size: #V9027-1 50ug #V9027-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Clostridium botulinum HA-33/Hemagglutinin 33 Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0038
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Clostridium botulinum
Conjugates	Unconjugated
Target Name	HA33, Hemagglutinin 33 protein, HA-33, Non-toxin haemagglutinin HA34, Clostridium botulinum
Accession No.	Q45871
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Clostridium botulinum
HA-33/Hemagglutinin 33 Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.