

Anti-Clostridium perfringens pfo/Perfringolysin O/Theta-toxin Recombinant Antibody

Catalog No: #V9030



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9030-1 50ug #V9030-2 100ug


Description

Product Name	Anti-Clostridium perfringens pfo/Perfringolysin O/Theta-toxin Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0014
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization, WB
Species Reactivity	Clostridium perfringens, Streptococcus dysgalactiae
Conjugates	Unconjugated
Target Name	Perfringolysin O, PFO, Theta-toxin, Thiol-activated cytolysin, pfo, pfoA, pfoR, CPE0163, Streptolysin O, SLO, Thiol-activated cytolysin, slo, Clostridium perfringens (strain 13 / Type A), Streptococcus dysgalactiae subsp. equisimilis (Streptococcus equisimilis)
Accession No.	P0C2E9, Q54114
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.


Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Clostridium perfringens pfo/Perfringolysin
O/Theta-toxin Antibody .



Western blot analysis was performed using anti-Theta-toxin monoclonal antibody at 1ug/mL on various samples. Lane 1: recombinant *Clostridium perfringens* pfo/Perfringolysin O/Theta-toxin Lane 2: negative control

Note: This product is for in vitro research use only and is not intended for use in humans or animals.