

Anti-Bacteriophage H19B stxB/SLT-1b/VT1 Recombinant Nanobody

Catalog No: #V9032



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9032-1 50ug #V9032-2 100ug


Description

Product Name	Anti-Bacteriophage H19B stxB/SLT-1b/VT1 Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	ST0056
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column.
Applications	ELISA, Neutralization, SPR
Species Reactivity	Enterobacteria phage H19B (Bacteriophage H19B)
Conjugates	Unconjugated
Target Name	Shiga-like toxin 1 subunit B, SLT-1 B subunit, SLT-1b, SLT-1b, Verocytotoxin 1 subunit B, Verotoxin 1 subunit B, stxB, sltB
Accession No.	P69179
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Bacteriophage H19B stxB/SLT-1b/VT1 Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.