

Anti-Escherichia coli Shiga toxin 2 subunit B/stx2eB
Recombinant Nanobody

Catalog No: #V9036



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9036-1 50ug #V9036-2 100ug


Description

Product Name	Anti-Escherichia coli Shiga toxin 2 subunit B/stx2eB Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	ST0036
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column.
Applications	ELISA
Species Reactivity	Escherichia coli
Conjugates	Unconjugated
Target Name	stxB2e, Shiga toxin 2e B-subunit, Stx2B, SLT-IIeB, Escherichia coli
Accession No.	Q47644
Formulation	0.94 mMB Citric acid monohydrate, 16.7mMB Trisodium citrate dihydrate, 6.2%B Sucrose, pHB 6.5, 0.01%B Tween 80.
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Escherichia coli Shiga toxin 2 subunit B/stx2eB Nanobody .

Note: This product is for in vitro research use only and is not intended for use in humans or animals.