

Anti-Staphylococcus aureus Alpha-toxin/hly/AT Recombinant Antibody

Catalog No: #V9039

Package Size: #V9039-1 50ug #V9039-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Staphylococcus aureus Alpha-toxin/hly/AT Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0055
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, Neutralization, SPR
Species Reactivity	Staphylococcus aureus
Conjugates	Unconjugated
Target Name	hly, hla, Alpha-hemolysin, Alpha-toxin, Alpha-HL, AT, Staphylococcus aureus, Staphylococcus aureus (strain NCTC 8325 / PS 47)
Accession No.	P09616, Q2G1X0
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.


Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Staphylococcus aureus Alpha-toxin/hly/AT Antibody



Western blot analysis was performed using anti-a1AT/Serpin A1/SERPINA1 monoclonal antibody at 1ug/mL on various samples. Lane 1: recombinant human a1AT/Serpin A1/SERPINA1 Lane 2: negative control

Note: This product is for in vitro research use only and is not intended for use in humans or animals.