

Anti-Microcystis aeruginosa Microcystin-LR Recombinant Antibody

Catalog No: #V9042

Package Size: #V9042-1 50ug #V9042-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Microcystis aeruginosa Microcystin-LR Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0001
Isotype	IgG2a
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA, WB
Species Reactivity	Microcystis aeruginosa
Conjugates	Unconjugated
Target Name	microcystin LR, Microcystin-LR, Microcystis aeruginosa toxin, toxin-LR, 1, 7-anhydro[D-alanyl-L-leucyl-(3S)-3-methyl-D-β-aspartyl-L-arginyl-(2S, 3S, 4E, 6E, 8S, 9S)-3-amino-4, 5, 6, 7-tetradehydro-9-methoxy-2, 6, 8-trimethyl-10-phenyldecanoyl-D-γ-glutamyl-2, 3-didehydro-N-methylalanine], cyanoginosin LR, cyclo[2, 3-didehydro-N-methylalanyl-D-alanyl-L-leucyl-erythro-3-methyl-D-β-aspartyl-L-arginyl-(2S, 3S, 4E, 6E, 8S, 9S)-3-amino-4, 5, 6, 7-tetradehydro-9-methoxy-2, 6, 8-trimethyl-10-phenyldecanoyl-D-γ-glutamyl], Microcystis aeruginosa
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Microcystis aeruginosa Microcystin-LR Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.