

Anti-Microcystis aeruginosa Microcystin-LR Recombinant Nanobody

Catalog No: #V9043

Package Size: #V9043-1 50ug #V9043-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Microcystis aeruginosa Microcystin-LR Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	ST0016
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column.
Applications	ELISA
Species Reactivity	Microcystis aeruginosa
Conjugates	Unconjugated
Target Name	microcystin LR, Microcystin-LR, Microcystis aeruginosa toxin, toxin-LR, 1, 7-anhydro[D-alanyl-L-leucyl-(3S)-3-methyl-D-β-aspartyl-L-arginyl-(2S, 3S, 4E, 6E, 8S, 9S)-3-amino-4, 5, 6, 7-tetradehydro-9-methoxy-2, 6, 8-trimethyl-10-phenyldecanoyl-D-γ-glutamyl-2, 3-didehydro-N-methylalanine], cyanoginosin LR, cyclo[2, 3-didehydro-N-methylalanyl-D-alanyl-L-leucyl-erythro-3-methyl-D-β-aspartyl-L-arginyl-(2S, 3S, 4E, 6E, 8S, 9S)-3-amino-4, 5, 6, 7-tetradehydro-9-methoxy-2, 6, 8-trimethyl-10-phenyldecanoyl-D-γ-glutamyl], Microcystis aeruginosa
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Microcystis aeruginosa Microcystin-LR Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.