

Anti-Alternaria species Tenuazonic acid Recombinant Nanobody

Catalog No: #V9045

Package Size: #V9045-1 50ug #V9045-2 100ug

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Alternaria species Tenuazonic acid Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	ST0030
Isotype	VHH-8His-Cys-tag
Purification	Purified by Nickel column.
Applications	ELISA
Species Reactivity	Alternaria
Conjugates	Unconjugated
Target Name	TeA, Tenuazonic acid, AAC-toxin, Alternaria mycotoxin, Alternaria species.
Accession No.	ChEBI: 9440
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Alternaria species Tenuazonic acid Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.