

# Anti-Alternaria species Tenuazonic acid Recombinant Nanobody

Catalog No: #V9047



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #V9047-1 50ug #V9047-2 100ug


## Description

|                    |                                                                                                                                                                           |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name       | Anti-Alternaria species Tenuazonic acid Recombinant Nanobody                                                                                                              |
| Host Species       | Alpaca                                                                                                                                                                    |
| Clonality          | Monoclonal                                                                                                                                                                |
| Clone No.          | ST0052                                                                                                                                                                    |
| Isotype            | VHH-8His-Cys-tag                                                                                                                                                          |
| Purification       | Purified by Nickel column.                                                                                                                                                |
| Applications       | ELISA                                                                                                                                                                     |
| Species Reactivity | Alternaria                                                                                                                                                                |
| Conjugates         | Unconjugated                                                                                                                                                              |
| Target Name        | TeA, Tenuazonic acid, AAC-toxin, Alternaria mycotoxin, Alternaria species.                                                                                                |
| Accession No.      | ChEBI: 9440                                                                                                                                                               |
| Formulation        | 0.01M PBS pH 7.4                                                                                                                                                          |
| Storage            | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt. |

## Application Details

Purity:>95% as determined by SDS-PAGE.

## Images



SDS-PAGE for Anti-Alternaria species Tenuazonic acid Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.