

Anti-Snake Venom Neutralizing Recombinant Antibody

Catalog No: #V9049



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9049-1 50ug #V9049-2 100ug


Description

Product Name	Anti-Snake Venom Neutralizing Recombinant Antibody
Host Species	Human
Clonality	Monoclonal
Clone No.	ST0015
Isotype	IgG1
Purification	Protein A/G purified from cell culture supernatant.
Applications	IP, SPR, in vivo neutralization
Species Reactivity	General
Conjugates	Unconjugated
Target Name	short-chain neurotoxin
Accession No.	P01416
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>90% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Snake Venom Neutralizing Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.