

Anti-Clostridium botulinum botA/BOTOX Recombinant Antibody

Catalog No: #V9052

Package Size: #V9052-1 50ug #V9052-2 100ug



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Anti-Clostridium botulinum botA/BOTOX Recombinant Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	ST0004
Isotype	VHH-mFc
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Clostridium botulinum
Conjugates	Unconjugated
Target Name	Botulinum neurotoxin type A, BoNT/A, Bontoxilysin-A, BOTOX, botA, bna
Accession No.	P0DPI0
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Clostridium botulinum botA/BoNT Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.