

Anti-Clostridium botulinum botA/BOTOX Recombinant Nanobody

Catalog No: #V9054



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #V9054-1 50ug #V9054-2 100ug


Description

Product Name	Anti-Clostridium botulinum botA/BOTOX Recombinant Nanobody
Host Species	Alpaca
Clonality	Monoclonal
Clone No.	ST0041
Isotype	VHH-His
Purification	Protein A/G purified from cell culture supernatant.
Applications	ELISA
Species Reactivity	Clostridium botulinum
Conjugates	Unconjugated
Target Name	Botulinum neurotoxin type A, BoNT/A, Bontoxilysin-A, BOTOX, botA, bna
Accession No.	P0DPI0
Formulation	0.01M PBS pH 7.4
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Application Details

Purity:>95% as determined by SDS-PAGE.

Images



SDS-PAGE for Anti-Clostridium botulinum botA/BoNT Nanobody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.